

Encapsulated Bacteria Session 3: Sequence-based typing

Genomics and Clinical Microbiology 2024

Martin Maiden, Made Krisna, Kasia Parfitt, Keith Jolley

Department of Biology



UNIVERSITY OF
OXFORD

Clinical Scenario Results

Specimen	Patient age	Presentation	Organism	Culture and CSF?	Serogroup
2	15 months	Meningitis	<i>N. meningitidis</i>	Yes	
5	3 years (46 months)	Meningitis	<i>N. meningitidis</i>	CSF only	

Clinical Scenario Results

Specimen	Patient age	Presentation	Organism	Culture and CSF?	Serogroup
2	15 months	Meningitis	<i>N. meningitidis</i>	Yes	B
5	3 years (46 months)	Meningitis	<i>N. meningitidis</i>	CSF only	B

Notes


- These two patients are related.
- They have both received Meningococcal serogroup C conjugate (MCC) vaccines, according to the Irish immunisation schedule at the time of birth.

Discussion points

- What actions are you going to take?
- What additional information do you need, and how are you going to get it?


Phenotypic characterisation

Gram's Stain



Anti-microbial
resistance
testing


Serology



Metabolic
phenotyping
For species
identification



MALDI-TOF phenotyping



Matrix-assisted laser desorption/ionization-time of flight (MALDI-TOF) mass spectrometry.


- Enables rapid and inexpensive (per sample) species identification
- Requires:
 - bacterial sample (usually an isolate);
 - the equipment;
 - a database linking spectra to bacterial species identification.
- Commercial systems are widely installed and used in clinical laboratories.

Hou, T. Y., Chiang-Ni, C. & Teng, S. H. (2019). Current status of MALDI-TOF mass spectrometry in clinical microbiology. *J Food Drug Anal.* **27**, 404-414.

Nucleotide sequencing and clinical microbiology

- **Definitive:**
 - fundamental level of information;
 - any part of the genome can be accessed.
- **Reproducible:**
 - nucleotide sequences are either right or wrong and can be checked;
 - reverse mutations are (usually) rare.
- **Scalable:**
 - nucleotide sequencing technology can be conducted on one or many samples and on a few base pairs or a whole genome.
- **Manipulability:**
 - nucleotide sequences can be analysed with model-based methods.
- **Can be done from a PCR reaction or microbiome sample:**
 - YOU DON'T (necessarily) NEED AN ISOLATE!


```
1 gagttttat c gct tccatg a cgc agaagt t aac actttc g gatattct g atg agtcga a
61 aaa ttattct t gat aaagca g gaattact c tgc ttgttt a cga attaaa t cga agtcga c
121 tgc tggcgg a aaa tagaaa a att cgacct a tct ttgccc a gct cgagaa g ctc ttactt t
181 gcg accctt c gcc atcaat c aac gattct g tca aaaa c acg cgttga t tga ggaaga g
241 tgg cttaat a tgc ttggca c gtt cgtcaa g gac ttgttt a gat atgagt c aca ttgtt t
301 cat ggtaga g att ctctgt t tga catttt a aaa gagcgt g gat tactat c tga gtcga t
361 gct gtccaa c cac taatag t taa gaattc a tga gtaagt t tac tgaaca c tcc gtacct t
421 tcc agacgc c ttt ggctct t att aagctc a ttc agcctt c tgc cgtttt g gat ttaacc g
481 aag atgatt t cga ttctct g acg agtaac a aag ttgtga t tgc tactga c cgc tctcgt g
541 ctc gtcgtc g cgt tgaggt t tgc gtttat g gta cgtcga c ctt ttggga t taa cctcgt t
601 ttc ctgctc c tgt tgagtt t att gctcgc c tca ttgctt a tta ttgtca c ccc gtaaca a
661 ttc aaacgg c ctg tctcat c atg gaagcg g ctg aattta c gga aaacat t att aatgac g
721 tgc agcgtc c ggt taaagc c gct gaattt t tgc cgttta c ctt cgtgtt a cgc gcagga a
781 aca ctgacg t tct tactga c gca gaagaa a acg tgcgtc a aaa attact g cgc gaagga g
841 tga tgaatg t tct aaaggt a aaa aacgtt c tgc cgtcgt c cct ggtcgt c cgc agcgtt t
901 gcg aggtac t aaa ggaagc c gta aagcgc c tgc tctttg t tat gtaggt g gta aacat t
961 tta attgca g ggg ctctcg c cct ttactt g agg ataat t atg tcta a ttc aaactg g
1021 cgc cgagcg t atg cgcgat c acc ttctcc a tct tggctt c ctt cgttgt c aga ttgtcg g
1081 tct tattac c att tcaact a ctc cgttat c cgc tggcga c tcc ttctga t tga agcgtt g
1141 tgg cgtctc c cgt cttctc c cat tgcgtc g tgg ccttgt t cag tggatt a agt tcatga a
1201 cag tggatt a agt tcatga a act gggtat a ttg accatg c
cct aagcat t tgt tcaag g
atg cctgac c gta ccaggg c
tgc cgttgt c gcc atctcaa a
tct cgcaca a tga cgaatt c
gct aatttg c ata ctgaca c
tct tcaatt g gag ttaaaa c
tct aatctc t ggc catctg a
ttt tctcgt c gtt tctaac a
cat ggaact a tgt ttactt t
cag tacttt a acg ataaag g
tat ggaact t tgc cgcgcg c
tct aagaag t tta agattg c
cct gcttat c acc ttcttg a
cga gaacgc g tac ttattc g
ttg cagtga a ata gtcagt t
cgc gattca a tca tgaatt c
taa aaattt t aat ttgttc c
tga aggttt t aat catggt t
ata agctgt t tct cacttc t
aag ctacat c gtc aacgtt a
ttc ttcaat g cat tcaagt g
gtg ctgata t tgc ttgtta t
att cttcgt t tcc gactac c
atg atggtg g tta ttatac c
cgg gcaata a cgt ttatgt t
cgc gattgt t ttc gctgaa t
gag gtgatt t atg ttgtgt g
tat gctaaa t tgc ttgtga g
gct tgcctc c gat aacaa t
agg ctctaa t gtt cctaac c
taa agctgt t aaa ggaatt c
taa gtgctt t gat ttgtgt g
tgc gattat t ctt gctggt g
tgc ttccct c gct ggtatg g
gca actgga c aat cagaaa g
tgg cattca g tgc ggcact t
gat gcttgt t tat caacga a
aaa tcttcc c aag caaaga c
aac ggcgtt t cag tatttt a
tga ggttga c tta gttcat c
cgc tactga a aag gatatt t
ttt tcatgt t att gataaa g
ggt gtaggt t att gctcta a
att gtagt t ttc atgctt c
ttc tgaatt t cag gctgat t
tga ggtgtt g ggc atttt a
gga taaccg c atc aagctt c
cga taactgt t gat atgatt g
tgt atctgt t act gagaag t
act tgaatt t aat aacact a
gaa cgagaa g acg gttacg c
tat tgcgta t gag tataat t
gct gggggc c tcc actatga a
aat ggcaga g gct catgt g
cga cagatt a gag gctttt t
tgg tctgat g gtt cttgct g
tca gatttt t tgt gtcgct g
gca ctttat g cga aacact c
caa tgcgcg c cag ttaaat a
cgc agttcg c tac acgcag g
taa aggtga g cgc cttaaa g
cgt taacaa a aag tcaagat a
caa ctcaat a aaa accaag c
cgc caactt c gga atgaaa a
act taacaa g ctg ggttac g
aag agagat g aga ttgag g
5281 tgg gaaaag t tac tttagc c gac gtttt g cgg cgcac c
5341 aat ttatgc g cgt tttagat a aaa atgatt g ggc tatcca a cct gca
```



Frederick Sanger (1918-2013)

Sanger, F., Air, G. M., Barrell, B. G., Brown, N. L., Coulson, A. R., Fiddes, C. A., Hutchison, C. A., Slocombe, P. M. & Smith, M. (1977). Nucleotide sequence of bacteriophage phi X174 DNA. *Nature* 265, 687-695.

Molecular approaches to clinical and public health microbiology



The role of genetics and genomics


- Nucleotide sequencing facilitates specimen characterisation,
 - sequence-based typing, up to the level of whole genome sequences (WGSs),
 - amplification techniques increase the amount of material for characterisation;
- Sequences are definitive, reproducible, and comparable.
 - easily stored, transported, and manipulated with analysis algorithms;
- However, datasets are large and complex,
 - requiring interpretation for the practitioner.

Steps in the characterisation of meningococci

- **Diagnosis:**
 - Discrimination from other organisms e.g. *Streptococcus pneumoniae* and *Haemophilus influenzae*.
- **Characterisation:**
 - Differentiation (e.g. meningococcus from other *Neisseria*);
 - Serogroup (conjugate polysaccharide vaccines available against serogroups A, C, Y, W, but not B);
 - Protein vaccine antigens (e.g. PorA, BAST typing), protein-based vaccines have been developed as substitutes to serogroup B vaccines;
 - The discrimination of isolates is essential for outbreak investigation;
 - Monitoring local and global trends in disease & the identification of likely epidemic strains.


Within species diversity *Neisseria meningitidis*

- Meningococci are highly diverse antigenically and genetically,
 - this diversity is structured.
- 12 capsular serogroups,
 - 6 associated with invasive disease.
- Extensive evidence of HGT, but stable lineages are present,
 - these are associated with phenotypes, including invasive disease.



Ganesh, K., Allam, M., Wolter, N., Bratcher, H. B., Harrison, O. B., Lucidarme, J., Borrow, R., de Gouveia, L., Meiring, S., Birkhead, M., Maiden, M. C., von Gottberg, A. & du Plessis, M. (2017). Molecular characterization of invasive capsule null *Neisseria meningitidis* in South Africa. *BMC Microbiology* **17**, 40.

Case Study, assembling a 'Strain type' or 'fine type': *Neisseria meningitidis* characterisation




Component	Phenotypic	Genotypic
Capsule	Serogroup	<i>cps</i> region
OMPS	Serotype, Subtype, etc.	<i>porA</i> , <i>porB</i> , <i>fetA</i> , etc.
Housekeeping genes	MLEE	MLST
Ribosomes	MALDITOF	16s rRNA, rMLST

***Neisseria meningitidis* C: P1.21-15,16: F1-7: ST-10217 (cc10217)**

Jolley, K. A., Brehony, C. & Maiden, M. C. (2007). Molecular typing of meningococci: recommendations for target choice and nomenclature. *FEMS Microbiol Rev* **31**, 89-96.


Meningococcal genome characterisation



Parkhill, J., Achtman, M., *et al.* Spratt, B. G. & Barrell, B. G. (2000). Complete DNA sequence of a serogroup A strain of *Neisseria meningitidis* Z2491. *Nature* **404**, 502-506.

Bratcher, H. B., Corton, C., Jolley, K. A., Parkhill, J. & Maiden, M. C. (2014). A gene-by-gene population genomics platform: *de novo* assembly, annotation and genealogical analysis of 108 representative *Neisseria meningitidis* genomes. *BMC Genomics* **15**, 1138.

Typing targets in the meningococcal genome



Antigen genes:

- Capsule (*cps*), serogroup;
- OMPs (*porA*, *porB*, *fetA*), serotype, subtype fine type;
- BAST antigens.

Housekeeping genes:

- Ribosomal RNA genes (*rrna*) genus identification;
- Ribosome protein genes (*rps*), species and sub species characterisation (rMLST);
- MLST (metabolism), subspecies identification
- cgMLST, high-resolution typing.

Antibiotic Resistance genes.


rplF assay

- **Problem:** rapid cost effective determination of *Neisseria* species from 1000s of isolates obtained in the MenAfriCar surveys.
 - 16s rRNA sequencing too cumbersome and insufficient discrimination or resolution.
- **Solution:** identify a short gene fragment (~400bp in length) diagnostic for species.
- **Implementation:** Phylogeny of *Neisseria* species generated and compared to phylogenies of individual genes
 - A fragment of *rplF* gene was congruent with clusters in the phylogeny of all genes and used for the assay.


Bennett, J. S., Watkins, E. R., Jolley, K. A., Harrison, O. B. & Maiden, M. C. (2014). Identifying *Neisseria* species using the 50S ribosomal protein L6 (*rplF*) gene. *J Clin Microbiol* 52, 1375-1381.

Comparison of 16s rRNA sequencing and *rplF* sequencing

16s rRNA gene fragment phylogeny (single rRNA encoding gene)




rMLST phylogeny (53 ribosomal protein encoding genes)




Bennett, J. S., Watkins, E. R., Jolley, K. A., Harrison, O. B. & Maiden, M. C. (2014). Identifying *Neisseria* species using the 50S ribosomal protein L6 (*rplF*) gene. *J Clin Microbiol* 52, 1375-1381.

Rapid species assignment: *rplF* sequence

rMLST phylogeny (53 genes)




rplF fragment (413bp) phylogeny





Bennett, J. S., Watkins, E. R., Jolley, K. A., Harrison, O. B. & Maiden, M. C. (2014). Identifying *Neisseria* species using the 50S ribosomal protein L6 (*rplF*) gene. *J Clin Microbiol* 52, 1375-1381.

Neisseria species distribution




Diallo, K., *et al.*, Greenwood, B. M. & Maiden, M. C. (2016). Pharyngeal carriage of *Neisseria* species in the African meningitis belt. *J Infect* **72**, 667-677.

Meningococcal genogrouping: a simpler problem



Harrison, O. B., Claus, H., Jiang, Y., Bennett, J. S., Bratcher, H. B., Jolley, K. A., Corton, C., Care, R., Poolman, J. T., Zollinger, W. D., Frasch, C. E., Stephens, D. S., Feavers, I., Frosch, M., Parkhill, J., Vogel, U., Quail, M. A., Bentley, S. D. & Maiden, M. C. J. (2013). Description and nomenclature of *Neisseria meningitidis* capsule locus. *Emerging Infectious Diseases* 19, 566-573.


Genogrouping: a PCR will do it



PorA sequence variation and serosubtyping

- Major outer membrane protein (OMP);
- Serosubtyping antigen (originally with mAbs) and a major vaccine component;
- Most antigenic variability in two variable regions (VRs) of the protein and gene.
- A single PCR amplification and sequencing reaction can provide high discrimination among isolates.


Structure of PorA P1.2 variant




Tzitzilonis, C., Prince, S. M., Collins, R. F., Achtman, M., Feavers, I. M., Maiden, M. C. & Derrick, J. P. (2006).
Structural variation and immune recognition of the P1.2 subtype meningococcal antigen. *Proteins* 62, 947-955.

Meningococcal characterization: sequence-based subtyping


- The two most variable parts of the protein (VR1 and VR2) are close together,
 - one sequencing reaction can determine the parts of the gene encoding both of these.
- *porA* gene fragment sequencing therefore extracts a lot of information for the amount of sequencing performed.



Sequencing *porA* VR1 and VR2



Direct sequencing: steps, reagents, equipment



Exercise

Specimen	Patient age	Presentation	Organism	Serogroup	PorA
2	15 months	Meningitis	<i>N. meningitidis</i>	B	
5	3 years (46 months)	Meningitis	<i>N. meningitidis</i>	B	

The laboratory has been able to obtain sequences of the *porA* gene (encoding the PorA protein) from each of the specimens:

- Specimen 2 - from DNA extracted from the culture;
- Specimen 5 – from direct sequencing of an amplicon of the *porA* gene obtained directly from the CSF.

Query the [PubMLST.org/neisseria](https://pubmlst.org/neisseria) website with the sequence data provided, to identify the PorA variants present in each of these meningococci.

Exercise -Results

Specimen	Patient age	Presentation	Organism	Serogroup	PorA
2	15 months	Meningitis	<i>N. meningitidis</i>	B	P1.7-2,4
5	3 years (46 months)	Meningitis	<i>N. meningitidis</i>	B	P1.7-2,4

Use the sequence data provided, along with the PubMLST.org/neisseria website, to identify the PorA variants present in each of these meningococci.

- Does this information reassure you, or not?
- Why?
- What should you do now (if anything)?

Encapsulated Bacteria: Web-based sequence querying

Genomics and Clinical Microbiology 2024

Martin Maiden, Made Krisna, Kasia Parfitt, Keith Jolley,
Department of Biology



UNIVERSITY OF
OXFORD

The plot thickens!

A further case of meningitis another family member has been recorded; however, this CSF sample yielded no further information.

A review of surveillance data shows four other cases of meningococcal disease occurred in this extended family over a two-year period:

Specimen 6 – Derived from CSF of 6 month-old with invasive meningococcal disease. This was confirmed as containing serogroup B meningococcal DNA by real-time PCR. There is no isolate for this specimen.

Specimen 7 – Derived from CSF from a 5 month-old child. This has been confirmed as containing serogroup B meningococcal DNA by Real-time PCR. An isolate has also been obtained from this specimen and stored.

Specimen 8 – Derived from CSF from a 5 month-old child. This has been confirmed as containing serogroup B meningococcal DNA by Real-time PCR. There was no isolate with this specimen.

Specimen 9 – Derived from CSF from a 10 month-old child. This has been confirmed as containing serogroup B meningococcal DNA by Real-time PCR. There was no isolate with this specimen.

Questions


- What is your reaction to this series of cases?
- What action may be required?
- What additional information to your need?
- What additional support might you need?

Multi locus typing approaches

While single locus typing (e.g. 16s rRNA sequencing) is useful in many circumstances there are advantages in typing multiple loci, including:

- Increased resolution;
- Additional information,
 - antimicrobial resistance, vaccine antigens, virulence determinants;
- Mitigating the impacts of recombination on isolate characterization (more on this later);
- Identification of microbial lineages.

Neisseria meningitidis characterisation




Component	Phenotypic	Genotypic
Capsule	Serogroup	<i>cps</i> region
OMPS	Serotype, Subtype, etc.	<i>porA</i> , <i>porB</i> , <i>fetA</i> , etc.
Housekeeping genes	MLEE	MLST
Ribosomes	MALDITOF	16s rRNA, rMLST

***Neisseria meningitidis* B: P1.7,16: F3-3: ST-32 (cc32)**

Jolley, K. A., Brehony, C. & Maiden, M. C. (2007). Molecular typing of meningococci: recommendations for target choice and nomenclature. *FEMS Microbiol Rev* **31**, 89-96.

Typing targets in the meningococcal genome



Antigen genes:



- Capsule (*cps*), serogroup;
- OMPs (*porA*, *porB*, *fetA*), serotype, subtype fine type;
- BAST antigens.

Housekeeping genes:

- Ribosomal RNA genes (*rrna*) genus identification;
- Ribosome protein genes (*rps*), species and sub species characterisation (rMLST);
- MLST (metabolism), subspecies identification
- cgMLST, high-resolution typing.

Antibiotic Resistance genes.

First generation genomic typing: single locus and MLST



B: P1.7,16: F5-1: ST-33 (cc32)


Antigen type
(fine type)

Sequence type &
clonal complex

Maiden, MCJ, Bygraves, JA, Feil, E, Morelli, G, Russell, JE, Urwin, R, Zhang, Q, Zhou, J, Zurth, K, Caugant, DA, Feavers, IM, Achtman, M & Spratt, BG. (1998). Multilocus sequence typing: a portable approach to the identification of clones within populations of pathogenic microorganisms. *Proc Natl Acad Sci USA* 95, 3140-3145.

MLST: allele-based analyses

- *Neisseria* seven-locus ST summarises 3,284bp.
- That is 0.15% of the 2.18Mbp genome.
- 14,199 STs in PubMLST database (20th January 2019).
- 702-1036 alleles per locus.
- Many polymorphisms per locus.



Russell, J. E., Urwin, R., Gray, S. J., Fox, A. J., Feavers, I. M. & Maiden, M. C. (2008). Molecular epidemiology of meningococcal disease in England and Wales 1975-1995, before the introduction of serogroup C conjugate vaccines. *Microbiology* 154, 1170-1177.

Exercise

- Use the sequence information supplied with the <https://pubmlst.org/neisseria> website to complete the following information for each specimen as far as possible:
 - MLST;
 - PorA antigen;
 - FetA antigen.
- Combine with the epidemiological and clinical information.
- What are your conclusions?

Exercise

[illegible]

Exercise

[illegible]

Extended meningococcal disease outbreak in Éire

specimen	Age (mo)	town	serogroup	<i>abcZ</i>	<i>adk</i>	<i>aroE</i>	<i>fumC</i>	<i>gdh</i>	<i>pdhC</i>	<i>pgm</i>	ST	cc	PorA	FetA
A (03/2010)	6	A	B	3	6		24	11	6			41/44	P1.7-2,4	F1-21
B * (11/2010)	5	A	B	3	6	9	24	11	6	9	6697	41/44	P1.7-2,4	F1-21
C (11/2010)	5	A	B	3	6	9	24	11	6	9	6697	41/44	P1.7-2,4	F5-12
D (01/2012)	10	B	B	4	5	2				20				
E (03/2013)	15	A	B	3	6	9	24	11	6	9	6697	41/44	P1.7-2,4	F5-12
F (04/13)	46	A	B		6	9	24	-	6			41/44	P1.7-2,4	
G (06/2013)	9	B	-	-	-	-	-	-	-	-	-	-	-	-
H* (11/2013)	6	A	B	3	6	9	24	11	6	9	6697	41/44	P1.7-2,4	F5-12

Discussion points

- What do you now understand about the relationships among these specimens?
- What additional information/actions are possible with these molecular typing data?
- What public health action, if any, is required?

How Sanger (dideoxy chain termination) sequencing works


Martin Maiden



UNIVERSITY OF
OXFORD


The remarkable contributions of Fred Sanger

- Linear structure of biological information.
- Protein sequencing,
 - insulin sequence.
- DNA sequencing,
 - dideoxy method.
- An efficient means of determining *phenotype* is inferring it from *genotype*:
 - and that means DNA sequence;
 - protein sequences can be deduced from DNA sequence.




Fred Sanger 1918-2013: Photo MRC

Di-deoxy chain-terminating nucleotides



Key components of Sanger sequencing (sequencing by synthesis)



The sequence ladder generated by cycle sequencing

3' -TACCGATAGCTGCTTTTG-5'

ATGGCt
ATGGCTa
ATGGCTAt
ATGGCTATc
ATGGCTATCg
ATGGCTATCGa
ATGGCTATCGAc
ATGGCTATCGACg

Deoxynucleotides:

A, C, G, T.

Dideoxynucleotides:

a, c, g, t.

Fluorescent Labels:



Detection of sequence ladders


- Separation of reaction products:
 - Gel based, e.g. ABI 377;
 - Capillary, e.g. ABI 3730.

This contrasts with ...

- Progressive detection of products:
 - e.g. Illumina platform.


Gel-based detection (...now essentially obsolete)

- Sequence reactions loaded on to an acrylamide slab gel.
- Separated by electrophoresis;
- Laser/camera system detects the products as they pass the end of the gel.



Capillary sequencing (still used)

- Sequence reactions loaded into a capillary filled with acrylamide.
- Separated by electrophoresis through the capillary
- Laser/camera system detects the products as they pass a window at the end of the capillary.
- Capillaries are automatically flushed-out and re-filled with polymer between samples.




Separation of cycle sequencing products



- Performed using automated sequencers
 - e.g. Applied Biosystems capillary sequencers
- Different sequencers have different capacities:
 - ABI 310 single capillary;
 - ABI 3100 16 capillary;
 - ABI 3070 48 or 96 capillary.


ABI 3730



Inside the ABI 3730



Visualisation of sequence traces



Sequencing by synthesis progresses one way ... and sometimes goes wrong

3' - TACCGATAGCTGCTTTTG - 5'

ATGGCt
ATGGCTt
ATGGCTAt
ATGGCTATc
ATGGCTATCg
ATGGCTATCGa
ATGGCTATCGAg
ATGGCTATCGACg

Deoxynucleotides:
A, C, G, T.

Dideoxynucleotides:
a, c, g, t.

Fluorescent Labels:


Errors can be mitigated by redundancy

- In a typical Sanger sequencing reaction there are hundreds (or thousands) of templates.
- These can be primed many times, increasing coverage.


Note:

- As the number of templates goes down, the possibility of random (non-systematic) error increases.
- Some errors are systematic – a ‘wrong’ base is likely to be incorporated in *in vitro* sequencing or the sequence is hard to read.

Accuracy by redundancy

Read 1	ATGGCTATCGACGAAAAC
Read 2	AT C GCTATCGACGAAAAC
Read 3	ATGGC A ATCGACGAAAAC
Read 4	ATGGCTATCGA A GAAAAC
Read 5	ATGGCTATCGACGAA A T C
Read 6	T TGGCTATCGACGAAAAC
Read 7	ATGGCTATC C ACGAA-AC
Read 5	ATGGCTATCGACGAAA A G
Read 6	ATGGCTA A CGACGAAAAC
Read 7	ATGGCTATCGACG T AAAC
Consensus	ATGGCTATCGACGAAAAC

Improved accuracy is achieved with sequencing both strands



The data from the reverse strand is **reverse complemented** (usually with a computer these days!) and used to generate a **consensus** sequence.

Factors affecting sequence accuracy

- Quality of the sequencing method.
- Number of templates in the reaction.
- Number of repeats.
- Strand coverage.

Note that although the best possible accuracy should be the goal, even poor quality data may be helpful clinically.

The quality of the data must be sufficient for the clinical need.