
Encapsulated Bacteria Session 2: Bacterial identification, Linux, and BLAST

Genomics and Clinical Microbiology 2024
Keith Jolley, Made Krisna, Kasia Parfitt, Martin Maiden,
Department of Biology



16S Ribosomal RNA sequencing


- Carl Woese initiated a unique research programme to obtain a tree of life.
- Ribosomal RNA was an accessible and tractable nucleic acid, which he characterised with thin layer chromatography sequencing.
- To his surprise, these analyses resulted in the discovery of a new domain of life first called the Archaebacteria (1977) and subsequently the Archaea (1990).



Sapp, J. & Fox, G. E. (2013). The singular quest for a universal tree of life. *Microbiol Mol Biol Rev.* **77**, 541-550.

16S rRNA gene analysis

- 16S rRNA an essential structural component of the ribosome,
 - highly conserved.
- Varies among organisms,
 - most variation is at the genus level.
- 16S rRNA gene, ~1.5Kbp, with:
 - universal, conserved regions, which can be used to generate primers;
 - nine variable regions which can be used as signatures to determine genus/species.



Yarza, P., Yilmaz, P., Pruesse, E., Glockner, F. O., Ludwig, W., Schleifer, K. H., Whitman, W. B., Euzeby, J., Amann, R. & Rossello-Mora, R. (2014). Uniting the classification of cultured and uncultured bacteria and archaea using 16S rRNA gene sequences. *Nat Rev Microbiol.* 12, 635-645.

16S rRNA sequencing in the clinic

Church, D. L., Cerutti, L., Gurtler, A., Griener, T., Zelazny, A. & Emler, **S.** (2020). Performance and Application of 16S rRNA Gene Cycle Sequencing for Routine Identification of Bacteria in the Clinical Microbiology Laboratory. *Clin* Microbiol Rev. 33, e00053-1

Sample Extraction

- · Verify method from various clinical samples/isolates being tested
- Use manual DNA extraction for up to 8-12 samples but automated method for high sample number
- Check DNA purity and concentration (i.e., NanoDrop)



FAST PCR

- Use High Tm (~60°C) 16S UMD primers + thermostable DNA polymerase
- Use FAST reagents used in FAST cycle-sequencing protocol (i.e., 1 h vs. 2.5 h for conventional PCR)
- Use FAST-ramping thermocycler



PCR Product Clean-Up


- ExoSAP-IT or equivalent enzymatic PCR product clean-up preferred
- Clean-up reaction removes unincorporated primers
- · Unincorporated nucleotides also degraded
- PCR product ready to use for sequencing without additional purification step (i.e., column purification kits)



FAST Cycle Sequencing

- Procedure similar to FAST PCR conditions
- Use BOTH forward (F) and reverse (R) primers = bidirectional read of rDNA
- · Reaction mix contains dNTPs and dye-labelled ddNTPs and decreased concentrations so extension is not inhibited
- Big Dye Xterminator extension product clean-up

16S rRNA 'species' id: process overview



Exercise: searching an NCBI database with 16SrRNA sequences

16s rRNA BLAST search online

- Navigate to NCBI BLAST home page:
 - https://blast.ncbi.nlm.nih.gov/Blast.cgi
- Click on:










- https://tinyurl.com/228hpuev
- Under 'Choose Search Set' click on 'rRNA/ITS databases'
- Check box 'Show results in a new window'. Click on









What bacterium is it?

- Look at the results of the search (scroll down the page as necessary).
- What is the likely organism that the 16S rRNA sequence came from?

Your Results

Sample	Candidate Bacterium
Test 1	Escherichia/Shigella???
Test 2	Mycobacterium tuberculosis
Test 3	Vibrio cholerae
Test 4	Neisseria mucosa
Test 5	Neissera gonorrheoae
Test 6	Neisseria meningitidis

Source sequences

Sample	Candidate Bacterium
Test 1	Escherichia coli (ENA MN900682)
Test 2	Mycobacterium tuberculosis (ENA AJ536031)
Test 3	Vibrio cholerae (ENA U10955)
Test 4	Neisseria lactamica (ENA FN995097.1)
Test 5	Neisseria gonorrheoae (ENA X07714)
Test 6	Neisseria meningitidis (ENA AJ239309)

Reflections

- How easy is it to get results?
- How easy is it to interpret them?
- How would you report them?
- How confident are you about the results?
- What information do you need to improve your confidence?

Practical: Introduction to Linux & BLAST