




UK Health  
Security  
Agency

# Short Read Sequencing Practical - 2025

Dr Jon Hubb  
Head of Colindale Sequencing Laboratory – Genomics Service Delivery and Quality Division


# Next Generation Sequencing methodologies

- Next Generation Sequencing
  - Catch all term
  - High throughput sequencing
  - Massively parallel sequencing
  - Terabytes of data
- Differing sequencing chemistries and engineering principles
- Sequencing of spatially separated, clonally amplified DNA.



# Patient sample pathways

- HCAI practical sessions:
  - Illumina short read sequencing (today)
  - Oxford Nanopore technologies long read sequencing (tomorrow)



# Preparation of Genomic DNA

Genomic DNA

Library preparation

Sequencing

- Purification of nucleic acid (genomic DNA and extra chromosomal)
- Extraction methods vary
  - One solution may not fit all
  - Often require pre-extraction preparation
- To name a few suppliers:
  - Qiagen - Qiasymphony or EZ1
  - Thermo - Kingfisher
  - Perkin Elmer - Chemagic
  - Roche - MagnaPure
- Which platforms is best?
  - Batch size
  - Sensitivity
  - Pre-extraction preparation
  - Cost
  - Run time
  - Spatial requirements
  - Established gold standard



# Library preparation


# Library Preparation (short read sequencing)

Genomic DNA

Library preparation

Sequencing


- Preparing the nucleic acid for sequencing
- DNA or RNA starting material
  - RNA needs to be reverse transcribed in cDNA
- Library preparation process
  - Fragmentation
  - End repair
  - dA-tailing 3'overhang
  - Ligation of adapter
  - PCR enrichment
  - Clean-up
  - Normalise
- Many commercially available kits for this process



# Illumina Nextera XT Method

- **Tagmentation (Section 2.0)**


- Starting material: Diluted Genomic or PCR amplified DNA (0.2ng/uL)
- Modification of DNA using transposomes enzymes
- Fragments DNA through enzymatic cleavage with addition of short adapter sequences that correspond to the sequencing primer sites
- Mix transposome, buffer and DNA and incubate for 5mins at 55°C
  - Reaction stopped by adding stopping buffer



# Illumina Nextera XT Method

## PCR Amplification (Section 3.0)

- PCR amplification of library
  - Incorporates terminal modification library with the incorporation of additional sequencer
- PCR primers are tripartite sequence structure
  - Complimentary adapter sequence: correspond to sequence primer binding site
  - Index: Adds a unique genetic tag onto the library for each samples eg.
    - UDP0001
      - i7 CGCTCAGTTC
      - i5 TCGTGGAGCG
  - Capture sequences: Perform PCR o modified DNA to include indices and P5 and P7 sequences







# Illumina Nextera XT Method

## PCR purification (Section 4.0)


- Purification of dsDNA using AMPure XP bead-based reagent
- Negatively charge DNA binds beads
- Purification of beads using magnet
- Wash impurities
- Low salt elution buffer dissociates DNA




# Illumina Nextera XT Method

## Normalisation (Section 5.0)

- Bringing all individual libraries in plate to same molarity prior to sequencing
- 4 main methods
  - Bead based:
    - Illumina library preparation (Nextera XT and Illumina DNA preparation)
  - Manual
    - Determine molarity using concentration of library and library size
  - Real time PCR
    - KAPA real time SYBR green method
    - Technically challenging
  - Enzymatic – Normalase
    - Enzymatic method
    - 2-step normalase procedure
    - Rapid, simple method



Normalisation method



Pool

# Nextera XT library preparation method



# SPRI beads



# Illumina Sequencing

# NGS sequencing technologies

Genomic DNA

Library preparation

Sequencing

454



GS Jr. GS FLX+

Amplification Method	Emulsion PCR on beads	
Chemistry	Synthesis (pyrosequencing)	
Read length (bp)	400	700
Yield/run (Gb)	0.05	0.9
Primary Error	Indel	
Error rate	~1%	
Run time (hours)	10	20
Virus-related Publications	187	
Advantage(s)	Long reads, maturity	
Disadvantage(s)	Homopolymer misreads, high cost/Mb	

Illumina



MiSeq HiSeq

Amplification Method	Bridge PCR in situ	
Chemistry	Synthesis (reversible termination)	
Read length (bp)	250	125
Yield/run (Gb)	8	1,000
Primary Error	Substitution	
Error rate	~0.1%	
Run Time (hours)	39	276
Virus-related Publications	129	
Advantage(s)	Easy work flow, maturity	
Disadvantage(s)	Shortest reads, long run	

Ion  
Torrent



PGM Proton

Amplification Method	Emulsion PCR on beads	
Chemistry	Synthesis (H <sup>+</sup> detection)	
Read length (bp)	400	200
Yield/run (Gb)	2	10
Primary Error	Indel	
Error rate	~1%	
Run time (hours)	7	4
Virus-related Publications	13	
Advantage(s)	Low cost, fast run	
Disadvantage(s)	Homopolymer misreads	

PacBio



RS II

Amplification Method	No PCR	
Chemistry	Single-molecule real-time sequencing	
Read length (bp)	8,500	
Yield/run (Gb)	0.15	
Primary Error	Indel	
Error rate	~13%	
Run time (hours)	2	
Virus-related Publications	6	
Advantage(s)	Longest reads	
Disadvantage(s)	High error rate, expensive	

# Illumina Sequencers

Genomic DNA

Library preparation

Sequencing

- MiniSeq, MiSeq, NextSeq (550/1000), HiSeq 2500, NovaSeq 6000
- Related methods of sequencing
  - Subtle difference in clustering and sequencing chemistries
  - Vary in amount of data and run time considerably







# Which sequencer or kit?

Genomic DNA

Library preparation

Sequencing

- Overwhelming amount of Illumina sequencers and kits available
- Key questions:
  - What are you trying to do?
  - How big is your organism?
  - What coverage is required?
  - What's my budget?

	 iSeq 100*	 MiniSeq*	 MiSeq*†	 NextSeq 1000 & 2000*
Output Range	144 Mb - 1.2 Gb	1.65–7.5 Gb	0.3–15 Gb	30–360 Gb***
Run Time	9–19 hours	4–24 hours	5–55 hr	11–48 hours
Reads Per Run	4 million	7–25 million	1–25 million	100 million–1.2 billion***
Max Read Length	2 × 150 bp	2 × 150 bp	2 × 300 bp	2 × 150 bp
Samples Per Run§	1–8	50	1–384	30–90
Relative Price Per Sample§	Higher Cost	Mid Cost	Mid Cost	Low Cost
Instrument Price	Lowest Cost	Low Cost	Low Cost	Mid Cost

<https://emea.illumina.com/systems/sequencing-platforms/comparison-tool.html#/research-use-only/microbiology/small-whole-genome-sequencing>  
[https://emea.support.illumina.com/downloads/sequencing\\_coverage\\_calculator.html](https://emea.support.illumina.com/downloads/sequencing_coverage_calculator.html)




# Cluster Generation

Genomic DNA

Library preparation

Sequencing


- Flow cell has P5 and P7 oligonucleotide lawn
  - Flow cell has single stranded complementary sequences
- P5 and P7 sequences on library termini bind to complementary sequences on flow cell
- Immobilised to discrete locations on flow cell



# Cluster generation


Single-stranded molecule flips over and forms a bridge by hybridizing to adjacent, complementary primer

Hybridized primer is extended by polymerases



Double-stranded bridge is denatured



Result:  
Two copies of covalently bound single-stranded templates



- Bridge amplification - Amplification of these single stranded fragments allows immobilised cluster generation
- Cluster generation can be random or more organised (patterned)
- Cluster generation - ~1000 clonal library copies per cluster

# Sequencing By Synthesis

- Sequencing by synthesis
- All 4 nucleotides are included
- Labelled with a different fluorophore
- Contains a block on 3' carbon
- Each base incorporation = 1 cycle



4-Channel Chemistry				
	A	G	T	C
Image 1	●			
Image 2		●		
Image 3			●	
Image 4				●
Result	A	G	T	C



2-Channel Chemistry				
	A	G	T	C
Image 1	●		●	
Image 2	●			●
Result	A	G	T	C

1-Channel Chemistry				
	A	G	T	C
Image 1	●		●	
Image 2			●	●
Result	A	G	T	C

..... Intermediate chemistry step

# Sequence determination

- During sequencing emission of light to indicate which base has been incorporated
- Sequencer camera records the fluorescence emission from the flow cell
- Each cluster appears as pin-prick of light
- Cluster generation allows:
  - Amplification of library which improves signal acquisition
  - Reduces impact of phasing



# Illumina

